

## Introduction

The gene-editing of T cells has assisted in the creation of therapies such as Chimeric Antigen Receptor (CAR) T and T Cell Receptor (TCR) T-cell therapies. Using CRISPR-Cas9 ribonucleoprotein (RNP) complexes, these innovative therapies offer hope for patients with otherwise untreatable cancers and other devastating diseases. We have developed a microfluidic transfection device which allows volume exchange for convective transfection (VECT). The device uses non-viral means to efficiently deliver CRISPR-Cas9 payload through mechanical transfection, leading to high levels of TCR knock-out, while retaining normal T cell functionality.

## Objectives

- Establish operational ranges at which VECT can efficiently deliver *TRAC* CRISPR-Cas9 RNP complexes.
- Determine the effect of CRISPR-Cas9 delivery through VECT on cell functionality.
- Evaluate feasibility and efficiency of gene knockout when RNP components (gRNA and Cas9 nuclease), are delivered separately.
- Assess the efficiency of VECT mediated knockout when co-transfecting multiple RNP complexes.

## Experimental Design

Human primary T cells were isolated and activated 12 hours before transfection. Cells were washed and resuspended at a cell density of  $3 \times 10^6$  cells/mL. They were then combined with *TRAC* or *CD5* CRISPR-Cas9 RNP complexes at a concentration of  $18 \mu\text{g/mL}$  and transfected through microfluidic channels. Cells were plated in 1mL media and cultured for 120 hours. Gene-editing, and cell functionality was measured using flow cytometry.

## Results

The transfection of *TRAC* CRISPR/Cas9 RNP using VECT (Figure 1a.) resulted in a high editing efficiency (>80%) while maintaining high viabilities (up to 80%). We were also able to recover a majority of cells processed through the microfluidic system (80%). Our CRISPR/Cas9 mediated knockout retains a high editing efficiency at varying RNP concentrations (Figure 1b.  $6 \mu\text{g/mL}$  –  $30 \mu\text{g/mL}$ ) and cell densities (Figure 1c.  $1 \times 10^6$  –  $8 \times 10^6$  cells/mL).

Figure 1. Operational Range of *TRAC* CRISPR/Cas9 RNP mediated knockouts using VECT






Figure 2. Cell Functionality of TCR KO T cells generated with VECT



The VECT transfection system preserves cell functionality in *TRAC* knockout samples. T cell proliferation is not influenced by microfluidic transfection. Our microfluidic device has no impact on the immunophenotype of either the CD4+ or the CD8+ T cell population (Figure 2b.). The profile of four common genes associated with T cell exhaustion also remains unaffected by the VECT system (Figure 2c.).

Figure 3. Microfluidic serial processing and co-transfection of CRISPR/Cas9 RNPs in human T Cells

Figure 3a.  $5 \times 10^6$  T cells were processed with  $18 \mu\text{g/mL}$  of *TRAC* single guide RNA in the morning, plated in media and then transfected again with  $18 \mu\text{g/mL}$  of cas9 protein. Samples in "Condition 2" were processed at a higher pressure than "Condition 1." Figure 3b.  $3 \times 10^6$  T cells were transfected with our microfluidic device with  $18 \mu\text{g/mL}$  of *CD5* and *TRAC* RNP complex at the same time.



The VECT system has the ability to gene-edit T cells at a high efficiency when the components of the RNP complex are delivered separately (Figure 3a.). To test this capability, we delivered the *TRAC* gRNA in the morning through a microfluidic device, and then delivered the cas9 nuclease in the afternoon. Our transfection rates remained comparable to the RNP control where both components were delivered as a pre-complexed RNP payload. Additionally, when processing the *CD5* and *TRAC* CRISPR-Cas9 RNP complexes we saw a mean knockout of 25% co-transfection. (Figure 3b.).

## Conclusions

VECT delivery of the *TRAC* knockout shows a large operational range and high editing efficiency in human primary T cells. Our system preserves T cell functionality and can be used to both serially deliver and co-transfect RNP complexes.

## Future Plans

In the future we plan on developing the knock-in capabilities of VECT by employing a range of DNA templates for homology-directed repair together with CRISPR/Cas9 RNPs.